

Magnum Bio

FILTERS & MEMBRANES



Magnum Bio

Cat #	Diameter (mm)	Pore Size (µm)	Membrane	Price (INR)
SYRINGE FILTER (STERILE)				
Nylon013G045ILP	13	0.45	Nylon(Pall membrane)	66
Nylon013G020ILP	13	0.2	Nylon(Pall membrane)	66
Nylon025G045ILP	25	0.45	Nylon(Pall membrane)	66
Nylon025G020ILP	25	0.2	Nylon(Pall membrane)	66
PES013G045ILP	13	0.45	PES(Pall membrane)	66
PES013G022ILP	13	0.22	PES(Pall membrane)	66
PES025G045ILP	25	0.45	PES(Pall membrane)	66
PES025G022ILP	25	0.22	PES(Pall membrane)	66
PVDF013G045IL	13	0.45	PVDF (hydrophilic)	66
PVDF013G022IL	13	0.22	PVDF (hydrophilic)	66
PVDF025G045IL	25	0.45	PVDF (hydrophilic)	66
PVDF025G022IL	25	0.22	PVDF (hydrophilic)	66
PVDF013G045ILM	13	0.45	PVDF (Millipore membrane)	66
PVDF013G022ILM	13	0.22	PVDF (Millipore membrane)	66
PVDF025G045ILM	25	0.45	PVDF (Millipore membrane)	72
PVDF025G022ILM	25	0.22	PVDF (Millipore membrane)	72
PTFE013G045IL	13	0.45	Hydrophilic PTFE	60
PTFE013G022IL	13	0.22	Hydrophilic PTFE	60
PTFE025G045IL	25	0.45	Hydrophilic PTFE	66
PTFE025G022IL	25	0.22	Hydrophilic PTFE	66
PTFE013G045OL	13	0.45	Hydrophobic PTFE	60
PTFE013G022OL	13	0.22	Hydrophobic PTFE	60
PTFE025G045OL	25	0.45	Hydrophobic PTFE	60
PTFE025G022OL	25	0.22	Hydrophobic PTFE	60



Product Type	Grade	Diameter (mm)	Packing	Price (INR)
FILTER PAPER				
		80*80cm	100pcs/pack	7668
		60*60cm	100pcs/case	3600
		7cm	100pcs/box	180
		9cm	100pcs/box	240
	101	11cm	100pcs/box	360
Qualitative	102	12.5cm	100pcs/box	420
Filter Paper		15cm	100pcs/box	540
	103	18cm	100pcs/box	720
		20cm	100pcs/box	960
		24cm	100pcs/box	1440
		25cm	100pcs/box	1440
		28cm	100pcs/box	1920
		30*30cm	100pcs/pack	4800
			1000pcs/case	
		7cm	100pcs/box	600
		9cm	100pcs/box	
	201	11cm	100pcs/box	960
Ashless	202	12.5cm	100pcs/box	
Filter Paper		15cm	100pcs/box	
	203	18cm	100pcs/box	1440

Cat #	Product Description	Qty.	Price (INR)
MAG02PCR	0.2 ml PCR thin wall tubes, 1000/bag	1 bag	4800
MAGI607	0.2 PCR rack, 96 well	1 rack	1740
MAG567-M	0.5-1.5 ml, 80 well Reversible Microtube Rack	1	1080
MAG577-M	0.5-1.5 ml, 96 well Reversible Microtube Rack	1	1680
MAG577-C	0.5-1.5 ml, 96 well Reversible, Natural Cover	1	840
MAG-1111	96-well-Reversible Assorted including Cover	1 set	2400
MAG-023	4-Ways Microtube & Centrifuge Tubes Racks	1	1920
MAGP001	Floating water bath rack, 10 rack/bag	1 bag	4800
MAG-4999	0.1-10 ul Clear tip, 1000/bag, DNase, RNase free, 10 bag/cs	10000/cs	3600
MAG-2346	1-200 ul Yellow tip, 1000/bag, DNase, RNase free, 10 bag/cs	10000/cs	3600
MAG-1150	50-1000 ul Blue tip, 1000/bag, DNase, RNase free, 10 bag/cs	10000/cs	4800
MAG-5100	0.1 to 10 ul, clear tips, DNAs/RNase free sterile, rack pack	96/rack	1440
MAG-2402	1 to 200 ul, yellow tips, DNAs/RNase free sterile, rack pack	96/rack	1464
MAG-1152B	100 to 1300 ul long tips, DNAs/RNAs free, sterile rack pack	96/rack	2040
MAG-5109	0.1-10 ul Extra-long Filter tip, 96/rack, DNase, RNase free, 10 rack/cs	10 rack/cs	2640
MAG-5106	1-200 ul Filter tip, 96/rack, DNase, RNase free, 10 rack/cs	10 rack/cs	2640
MAG-5107	50-1000 ul Filter tip, 96/rack, DNase, RNase free, 10 rack/cs	10 rack/cs	2880
MAG-2057	1.5ml microcentrifuge tube, large cap	500/bag	3600
MAG-2911	2.0ml microcentrifuge tube, large cap	500/bag	5040
MAG-4999	Rack for 0.1-10 ul Clear tip, 10 racks/box	10 rack/box	780
MAG-2346	Rack for 1-200 ul Yellow tip, 10 racks/box	10 rack/box	780
MAG-1150	Rack for 50-1000 ul Blue tip, 6 racks/box	6 rack/box	780
MAG4131	1.6ml, 96 deep well plates, square, Individual pack, sterile	1 case (96units)	1140
MAG4231	2.2ml, 96 deep well plates, square, Individual pack, sterile	1 case (96units)	1140
MAG4100	96 well Silicone film	10 psc/bag	13200



Cat #	Product Description	Qty.	Price (INR)
MAG-F001	Film for 96 deep well	96 film/bag	6600
MAG2702	1.5 ml, 100 Well, Cardboard, white, 13.4*13.4*4.7 cm.	100/case	540
MAG-23-2265	15 ml PP, conical, Printed Sterile, 50/rack, 10rack/cs	500/cs	4200
MAG-23-2263	50 ml PP, conical, Printed Sterile, 25/bag, 20bag/cs	500/cs	4800
MAG-200	200 UL tips, Dnase RNase Endotoxin free pack of 1000, white	1000/cs	324
MAG-200-Y	200 UL tips, Dnase RNase Endotoxin free pack of 1000, yellow	1000/cs	324
MAG-1000	1000 UL tips, Dnase RNase Endotoxin free pack of 500, white	500/cs	204
MAG-1000-Y	1000 UL tips, Dnase RNase Endotoxin free pack of 500, yellow	500/cs	204
MAG-10 -LR	10 UL tips, Low retention, Dnase RNase Endotoxin free pack of 1000, white	1000/cs	720
MAG-200-LR	200 UL tips, Low retention, Dnase RNase Endotoxin free pack of 1000, white	1000/cs	600
MAG-1000-LR	1000 UL tips, Low retention, Dnase RNase Endotoxin free pack of 1000, white	500/CS	480

Range of Biological Stains & Dyes

Features:

- Ideally suited to molecular biology or microbiology applications
- Available in a wide range of innovative chemical packaging options.

Staining is an important technique used in microscopy to enhance contrast in the microscopic image. Stains and dyes are frequently used in biology and medicine to highlight structures in biological tissues.

The comprehensive range of Stains and dyes, which are frequently used in Microbiology, Hematology, Histology, Cytology, Protein and DNA Staining after Electrophoresis and Fluorescence Microscopy etc. Many of the stains and dyes have specifications complying certified grade of Biological Stain Commission, and suitable for biological research. Stringent testing on all batches is performed to ensure consistency and satisfy necessary specification particularly in challenging work such as histology and molecular biology.

Stains and dyes offered includes Dry – powder form Stains and dyes as well as wet - ready to use solutions.

Product Name	CAS No
ACRIDINE ORANGE	10127-02-3
ACRIFLAVIN	8063-24-9
ALCIAN BLUE 8GX	33864-99-2
ALIZARINE AR	72-48-0
ALIZARINE CYANINE GREEN	4403-90-1
AMARANTH	915-67-3
AMIDO BLACK 10B	1064-48-8
ANILINE BLUE (SPIRIT SOLUBLE)	28631-66-5
ANILINE BLUE (WATER SOLUBLE)	28983-56-4
AURAMINE	2465-27-2
AZOCARMINE G	25641-18-3
AZUR A	531-53-3
AZUR B	531-55-3
AZUR II	37247-10-2
AZUR II EOSINE	53092-85-6
BIEBRICH SCARLET	4196-99-0
BISMARCK BROWN®	5421-66-9
BISMARCK BROWN Y (G)	10114-58-6
BORAX CARMINE GRENACHER	-
BRIILLIANT CRESYL BLUE	81029-05-2
BRIILLIANT GREEN INDICATOR	633-03-4
BRIILLIANT YELLOW	3051-11-4
BROMOCRESOL PURPLE SODIUM SALT INDICATOR	62625-30-3
CARBOL FUCHSIN POWDER	4197-24-4
CARMINE	1390-65-4
CARMINIC ACID	1260-17-9
CARMOISINE A	3567-69-9
CELESTIN BLUE	1562-90-9
CHROMAZUROL S	1667-99-8

Product Name	CAS No
CHROMOTROPE 2R	4197-07-3
CHRYSOIDINE R	4438-16-8
CHRYSOIDINE Y	532-82-1
CONGO RED	573-58-0
COOMASIE BRILLIANT BLUE G 250	6104-59-2
COOMASIE BRILLIANT BLUE R 250	6104-58-1
CRESYL VIOLET (ACETATE)	10510-54-0
CRYSTAL VIOLET AR/ACS	548-62-9
CRYSTAL VIOLET	548-62-9
EOSIN BLUE	548-24-3
EOSIN SPIRIT SOLUBLE	6359-05-3
EOSIN YELLOW (WATER SOLUBLE)	17372-87-1
ERIOCHROME BLACK T	17372-87-1
ERIOCHROME BLACK T AR/ACS	17372-87-1
ERIOCHROME CYANINE R AR	3564-18-9
ERIOGLAUCINE	3844-45-9
ERITHROSIN B	16423-68-0
ETHYL VIOLET	2390-59-2
EVANS BLUE	314-13-6
FAST BLUE B SALT	14263-94-6
FAST GARNET GBC SALT Purified	101-89-3
FAST GREEN FCF	2353-45-9
FIELD STAIN A	-
FIELD STAIN B	-
FLUORESCEIN 95%	2321-07-5
FLUORESCEIN SODIUM	518-47-8
FUCHSIN ACID	3244-88-0
FUCHSIN BASIC	632-99-5
GIEMSA'S STAIN	51811-82-6
HEMATOXYLIN STAIN CERTIFIED	517-28-2

Range of Biological Stains & Dyes (contd...)

Product Name	CAS No
INDIAN INK	-
JANUS GREEN B	2869-83-2
JENNER'S STAIN	-
LEISHMAN'S STAIN	-
LIGHT GREEN (Light Green SF Yellowish)	5141-20-8
LISSAMINE GREEN	3087-16-9
M.T.T.TETRAZOLIUM AR	298-93-1
MALACHITE GREEN	2437-29-8
MARTIUS YELLOW	605-69-6
MAY AND GRUNWALD'S STAIN	-
METHYL BLUE	28983-56-4
METHYL GREEN	7114-03-6
METHYL VIOLET	8004-87-3
METHYLENE BLUE A	61-73-4
METHYLENE BLUE AR	61-73-4
METHYLENE BLUE POLYCHROME	-
NAPHTHOL GREEN B	19381-50-1
NEW FUCHSIN	3248-91-7
NIGROSINE ALCOHOL SOLUBLE	11099-03-9
NIGROSINE WATER SOLUBLE	8005-03-6
NILE BLUE CHLORIDE	2381-85-3
NILE BLUE SULPHATE	3625-57-8
NITRAZINE YELLOW	5423-07-4
OIL RED O	1320-06-5
ORACET BLUE 2R	4395-65-7
ORANGE G	1936-15-8
ORCEIN	1400-62-0
PARAROSANILINE BASE	467-62-9

Product Name	CAS No
PARAROSANILINE CHLORIDE	569-61-9
PATENT BLUE V (VF)	129-17-9
PHENOSAFRANINE	81-93-6
PHLOXIN B	18472-87-2
PONCEAU S	6226-79-5
PYRONIN G	92-32-0
PYRONINE B	2150-48-3
QUINOLINE YELLOW	8004-92-0
RHODAMINE B	81-88-9
ROSE BENGAL AR	632-69-9
RUTHENIUM RED	11103-72-3
SAFRANINE	477-73-6
SUDAN BLACK B	4197-25-5
SUDAN III	85-86-9
SUDAN IV	85-83-6
TARTRAZINE AR	1934-21-0
TETRAZOLIUM SALT AR	298-96-4
TETRAZOLIUM VIOLET AR	1719-71-7
THIOFLAVINE T Extra Pure	2390-54-7
TOLUIDINE BLUE	92-31-9
TROPAEOLIN O	547-57-9
TROPAEOLIN OO	554-73-4
TROPAEOLIN OOO	633-96-5
TRYPAN BLUE	72-57-1
VICTORIA BLUE	2580-56-5
WRIGHT'S STAIN	68988-92-1
XYLENE CYNOL FF	2650-17-1
XYLIDINE PONCEAU	3761-53-3