



# irisbiosciences

## Ordering Information

Request for separate catalogue

Ref	Contents	Amount
060201	<b>Minisys G</b> – 3D cell culture hydrogel	0.2 g / vial
060230	Human Stromal Vascular Fraction in cryopreservation medium	1x10 <sup>6</sup> cells
060231	Human Adipose mesenchymal Stromal/stem Cells in cryopreservation medium	1x10 <sup>6</sup> cells
060232	<b>Xeno-free</b> human Adipose mesenchymal Stromal/stem Cells	1x10 <sup>6</sup> cells
060233	<b>Xeno-free</b> human Stromal Vascular Fraction in xeno-free cryopreservation medium	1x10 <sup>6</sup> cells



060201



060233



060232

hSVF xeno-free cultured in 2D cell culture plasticware

hSVF xeno-free cultured in 3D minisys G

Osteogenesis  
Alizarin red staining of mineralized extracellular matrix (ECM)

Adipogenesis  
Oil red O staining of lipid vacuoles

Chondrogenesis  
Alcian blue staining of ECM proteoglycan (pellet cross-section)

